

# Polonium Analysis in Complex Samples

Andrew W. Nelson

*PhD/MPH Candidate*

*Schultz Group*


*Interdisciplinary Human Toxicology Program*

*University of Iowa, Iowa City, Iowa*


# Matrix of Hydraulic Fracturing Flowback Water/Produced Fluids

Analyte	Concentration (mg/L)
Chloride	146,667
Strontium	36,333
Sodium	29,333
Calcium	13,000
Barium	9,000
Magnesium	853
Manganese	3
Iron	43
Total Solids	277,666
Suspended Solids	783

# Polonium on SR Resin




Data from Eichrom Website



Data Provided by Dr. Dan McAlister 3

# Ca and Po Retention on TRU





TRU Resin

Data Provided by Dr. Dan McAlister

Analyte	Concentration (mg/L)
Chloride	146,667
Strontium	36,333
Sodium	29,333
Calcium	13,000
Barium	9,000

# Polonium on TRU



Data Provided by Dr. Dan McAlister

# Polonium on TRU HCl

- Can analyze complex samples
  - Flowback
  - Samples with >20,000 ppm Fe
- After TRU
  - Autodeposition
    - High recoveries >85%
  - CuS microprecipitation



# Suggestions to Eichrom

- Provide the Po data in HCl on the TRU website
- Add Po as one of the elements that can be analyzed by TRU

