

## **Disclaimer:**

**The presented information is intended solely for research purposes. It should not be seen or used as an endorsement for promoting the products been tested.**

**The Food and Drug Administration can neither encourage nor discourage the use of the products been presented.**

# **Observations from Analysis of Alpha/Beta Radionuclides in Food Using Quantulus 1220, Quantulus GCT, and Hidex 300 SL Liquid Scintillation Spectrometers**

Zhichao Lin, Stephanie Healey, Patrick Regan

**Winchester Engineering and Analytical Center  
Food and Drug Administration**

Denis Bonaccolta, Rick Troyer, Brad Ward,  
Richard Harazin, Michael Belobraydich

**Radiometric Solutions and Discovery Services  
PerkinElmer, Inc.**

# Presentation Outline

- Motivations
- Objectives
- Experimental Approach
- Sample Preparation
- Instrument Calibrations & Optimizations
- Observations & Discussions
- Summary

## ➤ Motivations

- Increasing challenges in safeguard FDA-regulated products from radioactive contamination call for a broad range of complementary radioanalytical techniques
- Detection and quantification of alpha and beta radioactivity are essential parts of FDA food safety compliance and emergency response programs
- Advancement of FDA radioanalytical capability and surge capacity hinges on comprehensive metrology knowledge
- Development of scientific collaborations and partnerships to promote and advance regulatory science

## ➤ Objectives

- To adapt the current FDA radioanalytical capability and surge capacity to more novel, efficient, and sensitive instrument technologies
- To identify the detection techniques that simplifies and expedites analysis of alpha and beta radionuclides in foods
- To ascertain performance characteristics of novel technology based liquid scintillation spectrometers
- To make effective use of instrumentations for different radioanalytical needs

## ➤ Experimental Approach

- Three liquid scintillation spectrometers different in detector designs and pulse analysis techniques were used
- The following performance characteristics were studied:
  - ✓ Background
  - ✓ Cerenkov Counting
  - ✓ Simultaneous  $\alpha/\beta$  Analysis
- The following  $\alpha/\beta$  radionuclides were analyzed:
  - ✓  $^3\text{H}$
  - ✓  $^{89}\text{Sr}$  and  $^{90}\text{Y}$
  - ✓  $^{90}\text{Sr}/^{241}\text{Am}$  and  $^{90}\text{Sr}/^{239}\text{Pu}$

# Detector Assembly & $\alpha/\beta$ Pulse Analyzer



- ✓ Low background Pb shield
- ✓ Pulse amplitude comparison (PAC)
- ✓  $\alpha/\beta$  discrimination based on pulse shape (PSA)
- ✓ Bismuth germanium oxide (BGO) detector guard
- ✓ Guard compensation technology (GCT) algorithm




- ✓ Low background Pb shield
- ✓ Guard detectors
- ✓ Pulse amplitude comparison (PAC)
- ✓  $\alpha/\beta$  discrimination based on pulse shape (PSA)



*Hidex 300 SL*


- ✓ Low background Pb shield
- ✓ Guard detector
- ✓ Triple PMT coincidence counting
- ✓  $\alpha/\beta$  discrimination based on pulse length index (PLI)
- ✓ Triple-to-Double Coincidence Ratio (TDCR)

# Quantulus 1220 LSC





# Quantulus GCT LSC




GCT Algorithm

Quantulus GCT Active & Passive shielding, BGO guard, & GCT algorithm for background reduction

# TDCR Based LSC – Hidex 300 SL

## Triple to Double Coincidence Ratio Counting (TDCR)



## ➤ Sample Preparation

- Standard Preparation

### **For Analysis of $^3\text{H}$ :**

10 mL  $\text{H}_2\text{O}$

10 mL of Ultima Gold LLT

Teflon-coated plastic vial

### **For Cerenkov Counting of $^{89}\text{Sr}$ , $^{90}\text{Y}$ :**

10 mL 0.1M HCl

Teflon-coated plastic vial

### **For Simultaneous Counting of $^{90}\text{Sr}/^{241}\text{Am}$ , $^{90}\text{Sr}/^{239}\text{Pu}$ :**

3 mL of 0.5M HCl

17 mL of Ultima Gold AB

Low  $^{40}\text{K}$  glass vial



## ➤ Instrument Calibrations & Optimizations

	<sup>3</sup> H Efficiency %	Background CPM	E <sup>2</sup> /B	<sup>89</sup> Sr Efficiency %	Background CPM	E <sup>2</sup> /B
<b>Quantulus GCT (<sup>3</sup>H window = 0 – 18.6 keV, <sup>89</sup>Sr window = 0 – 20 keV )</b>						
GCT=Off*	19.04	1.50	242	32.04	1.95	526
GCT=Low	22.40	0.68	738	na	na	na
GCT=High	22.43	0.50	1006	31.82	1.0	1013
<b>Quantulus 1220 (<sup>3</sup>H window = 5 – 320 channel, <sup>89</sup>Sr window = 20 – 350 channel)</b>						
	21.55	2.57	180	41.64	1.56	1111
<b>Hidex 300 SL (<sup>3</sup>H window = 5 – 250 channel, <sup>89</sup>Sr window = 5 – 350 channel)</b>						
TDCR	12.71	3.43	47	na	na	na
TDCR	13.06	6.56	26	na	na	na

\* All GCT=Off data shown uses Low Level Count Mode (LLCM)

## ➤ Observations & Discussions


A significant background reduction was observed in the spectral region most critical for  $^3\text{H}$  analysis.



Benefited by lower background, the measurement accuracy for low level  $^3\text{H}$  analysis was noticeably improved with Quantulus GCT when the GCT algorithm was applied.

# ➤ Observations on $^3\text{H}$ Analysis





## Quantulus GCT Water Results:



## Quantulus GCT Food Results:


Food Matrix	Known, Bq/g	Measured Bq/g
Spiked Whole Milk	$13.72 \pm 0.44$ (2s)	$13.42 \pm 0.51$ (2s)
Spiked Fat-Free Milk	$13.72 \pm 0.44$ (2s)	$13.62 \pm 0.50$ (2s)

# ➤ Observations on Alpha/Beta Spillover



*Quantulus GCT*

# Visual Presentation of $\alpha/\beta$ Puls Separation with Hidex 300 SL





## ➤ Observations on Alpha/Beta Spillover (con't)

LSC Counter	Alpha Counted as Beta, %	Beta Counted as Alpha, %	$\alpha/\beta$ Discriminator Setting
$^{90}\text{Sr}/^{239}\text{Pu}$			
Quantulus GCT	1.70	1.99	72
Quantulus 1220	0.65	0.57	52
Hidex 300 SL	3.73	5.34	3
$^{90}\text{Sr}/^{241}\text{Am}$			
Quantulus GCT	1.00	0.99	97
Quantulus 1220	0.34	0.38	56
Hidex 300 SL	2.85	5.34	3

$$^{90}\text{Sr } E_{\text{max}} = 0.546 \text{ MeV}$$

$$^{90}\text{Y } E_{\text{max}} = 2.28 \text{ MeV}$$

$$^{241}\text{Am } E_{\alpha} = 5.49 \text{ MeV}$$

$$^{239}\text{Pu } E_{\alpha} = 5.15 \text{ MeV}$$

## ➤ Observations on Cerenkov Counting $^{89}\text{Sr}$ , $^{90}\text{Y}$

Sample ID	Measured, Bq	Known, Bq	Difference, %
Quantulus GCT			
SR89A1	0.948	0.948	0.02
SR89B2	1.882	1.896	-0.73
SR89C3	2.847	2.844	0.12
SR89D4	3.844	3.791	1.39
SR89E1	0.082	0.095	-13.6
Quantulus 1220			
SR89A1	0.946	0.948	-0.15
SR89B2	1.918	1.896	1.18
SR89C3	2.828	2.844	-0.55
SR89D4	3.777	3.791	-0.37
SR89E1	0.107	0.095	12.4

## ➤ Summary

- The GCT algorithm works reliably for analyzing low-level  $^3\text{H}$  when Quantulus GCT spectrometer is calibrated and adapted to local environment.
- Quantulus 1220 still holds the benchmark  $\alpha/\beta$  separation in simultaneous  $\alpha/\beta$  analysis but its superiority in background reduction is surpassed by the GCT technology.
- Hidex TDCR-based LSC provides unique advantage over traditional LSC in sample quench determinations.
- Lack of an option for fine-tuning  $\alpha/\beta$  discriminator excludes user from using the optimum discriminator setting

## ➤ Summary (con't)

- Quantulus GCT LSC allows user to reprocess the spectrum that had already been collected with a new set of parameters without recounting the samples
- Hidex 300 SL uses Excel-based data file which can be easily used by user's Excel spreadsheet for automated sample activity calculation
- The application software for Quantulus GCT LSC was found to be easy to use. The entire well-known TriCab's menu layout is kept as it was with a few new add-on options
- Hidex 300 SL allow almost quantitative activity concentration determinations without standardization



## Acknowledgment

We would like to thank

Debanond Chakraborty  
Brian D. Scherer

Radiochemistry Unit  
New Hampshire Division of Public Health Services

For their assistances and use of Hidex 300 SL throughout this study



Questions?