

## SEPARATION AND ANALYSIS OF Sr-90 AND Zr-90 FOR NUCLEAR FORENSIC APPLICATIONS

Dominic Larivière, Radioecology Laboratory, Chemistry Department, Laval University, Québec City, Canada, dominic.lariviere@chm.ulaval.ca

Ana Paula Zattoni, Radioecology Laboratory, Chemistry Department, Laval University, Québec City, Canada.

In this work, a technological development to determine the age of radioactive strontium sources through the [Zr-90]/[Sr-90] ratio using mass spectrometry and liquid scintillation to quantify both isotopes is presented. Because Sr-90 and Zr-90 are isobaric interferences in mass spectrometry, a radiochemical separation to isolate Zr-90 has been shown to be mandatory prior to analysis. Four commercial resins (AG50W-X9, Dowex1-X8, Sr and DGA resins) were tested to isolate Zr-90 from Sr-90. Best performance was observed for the DGA resin, including recoveries higher than 99% for Zr-90. DGA has also demonstrated to be the faster approach and the most efficient not only to eliminate isobaric interferences from Sr-90, but also from Y-90, potentially present in samples containing high levels of radioactivity.